



Puritan Bennett™
840 Ventilator



COVIDIEN

positive results for life™



The Link Between Sedation and ICU Outcomes

Keeping patients comfortable and at ease while they're in the ICU can be very challenging. Patients often have limited consciousness or ability to communicate.^{1,2}


A full 71% of patients show signs of agitation at least once during their stay.¹ Out of compassion, clinicians often turn to sedation to relieve distress.¹


However, a growing body of research has confirmed a strong link between sedation and poor patient outcomes. When used inappropriately, sedation can lead to failure to wean, prolonged ICU stays and increased cost of care.¹

The Challenging Reality of Mechanical Ventilation

Although a necessary intervention, conventional modes[†] of mechanical ventilation are limited in their ability to properly manage a patient's work of breathing.³⁻⁵ In fact, in the ICU, 42% of increases in sedation are in response to patient-ventilator asynchrony.¹

Without a way to better manage work of breathing, increasing sedation can seem like the best option. But, increased sedation can prolong time on the ventilator.^{1,2}





In addition, conventional modes of ventilation are rigid in their delivery of a breath. This pattern can be at odds with the patient's natural breathing rhythm, facilitating the cycle of asynchrony, sedation and muscle weakness.^{3, 4}

Volume Controlled

If a patient wakes up and tries to participate in breathing, work of breathing (WOB) and/or anxiety can increase:

- WOB – Vt limit is not high enough or the breath is too short. This can also cause anxiety.
- Anxiety – Breath is too long
- Patient may show signs of agitation and is sedated³

Pressure Controlled

Relatively comfortable mode, patient can receive variable flow, but the patient does not have to do any work leading to an increased risk of muscle atrophy.^{7, 11}

Pressure Support

Spontaneous breathing is allowed, however, it's still a programmed breath. If the patient triggers the ventilator he or she receives the programmed breath.

The patient either has to work very hard or not at all if ventilator settings are not aligned.

- WOB – if the level of support is below patient demand, the WOB will increase
- Reduced activity – if the level of support exceeds patient demand, the patient may have periods of reduced activity leading to an increased risk of muscle weakness
- Ineffective triggers – found to be an independent predictor of longer mechanical ventilation duration,⁴ ineffective triggers can result in an increase in WOB and/or anxiety
- Patient may show signs of agitation and is sedated³

To the patient, waiting for that next breath feels like they're suffocating, causing panic. Being forced to breathe too quickly, too shallowly, or too deeply is unsatisfying and extremely tiring, leading to distress.⁵

Promote Natural Breathing

At Covidien, we believe mechanical ventilation can and should be more natural.


Our PAV™*+ software for the Puritan Bennett™ 840 ventilator is a breath type that better manages the patient’s work of breathing and promotes natural breathing compared to conventional mechanical ventilation†.8

PAV™*+ software manages the patient’s work of breathing differently than other current modes of mechanical ventilation† in the following ways.12

With PAV™*+ mode the patient defines rate, depth and timing.

- Flow is an indicator of demand. It tells us when the patient wants to begin inspiration, how deep the breath should be, when to end the breath and how often to breathe.
- PAV™*+ mode continuously measures patient demand by measuring flow and volume every 5 milliseconds.
- As patient demand changes, PAV™*+ mode changes support within the same breath.

Enabling the patient to define rate, depth and timing helps reduce excessive workload or patient agitation, potentially reducing the need for unnecessary sedation.13-17



PAV™*+ mode provides better synchrony with a patient’s breathing than pressure support ventilation modes (PSV)

Year	Author	Results
2011	Costa et al ¹³	“PAV+ improves patient-ventilator interaction, significantly reducing the incidence of end-expiratory asynchrony and increasing the time of synchrony.”
2009	Xirouchaki et al ¹⁴	“Compared to PS, PAV+ is associated with fewer interventions in terms of ventilator settings and sedative dose changes.”
2008	Xirouchaki et al ¹⁵	“Compared to PS, PAV+ increases the probability of remaining on spontaneous breathing, while it considerably reduces the incidence of patient-ventilator asynchronies.”
2007	Bosma et al ¹⁶	“PAV+ resulted in better quality of sleep compared to PSV.”



When the %Support is set, the patient and the ventilator are sharing the work of breathing as defined by the clinician.


- Work of breathing can be calculated using the equation of motion.¹⁸
- When R and E are known, it's possible to calculate patient-generated pressure (P_{musc}) and work of breathing in real time using the equation of motion.^{16, 18-21}

$$P_{MUSC} + P_{VENT} = (\text{flow} \times \text{resistance}) + (\text{volume} \times \text{elastance})$$

- PAV™*+ mode measures resistance and compliance every 4-10 breaths.
- Once %Support is set, clinicians can use the work of breathing (WOB) bar for real-time feedback on how much work the patient is doing.
- The work of breathing bar displays both total work of breathing (WOB_{TOT}) and the patient work of breathing (WOB_{PT}).
- Associated fatigue values for work of breathing are shown as being outside the green zone.

The work of breathing bar, when coupled with good clinical assessment, can help take the guesswork out of determining the appropriate level of mechanical ventilation support.

Providing real-time feedback on work of breathing enables the clinician to keep the patient at a sustainable level of work—reducing the risk for respiratory muscle atrophy, but off-loading enough work to avoid fatigue.^{7,9,11}



The work of breathing bar shows total work of breathing (WOB_{TOT}) and patient work of breathing (WOB_{PT})



VC		PA	V-TRIG	50 kg
V _T	V _{max}	% Supp	V _{max}	O ₂
365 ml	22 l/min	50 %	3.0 l/min	100 %
T _{res}	□	E _{max}	PEEP	
0.0	SQUARE	3 l/min	3.0 cmH ₂ O	

Clinician sets “%Supp” to define the patient’s WOB%age



A Good Fit for Your Institution

In addition to PAV™* software, the Puritan Bennett™ 840 ventilator features a full suite of software options, safety features and accessories to fit a variety of patients, from infant to adult.

Software Options

Leak Compensation: To help ensure patients are receiving the flow and volume they need, Leak Compensation software immediately detects and rapidly adjusts to changes caused by patient or airway interface movement, leaks and other barriers to targeted breath delivery. This helps prevent auto-triggering and as a result can help minimize patient-ventilator asynchrony.

In studies, the Puritan Bennett™ 840 ventilator with Leak Compensation software was shown to:^{5, 12}

- Synchronize to increasing and decreasing leaks in both obstructive and restrictive lung models and with PEEP 5 cm H₂O and 10 cm H₂O
- Require fewer breaths to synchronize, under all test conditions
- Perform equally in both invasive and noninvasive ventilation settings



NeoMode 2.0 Software: This software enables the clinician to adjust ideal body weight (IBW) without disconnecting the patient, thereby avoiding the additional risks associated with patient disconnection. The ventilator is also able to detect and compensate for patient leaks, reducing the number of nuisance alarms and improving patient safety.

Bi-Level Software: Permits spontaneous breathing at all times, reducing patient-ventilator asynchrony. This software supports biphasic or airway pressure release ventilation for extra flexibility. An active exhalation valve improves patient comfort and reduces ICU costs associated with continuous sedation.

Volume Control Plus: This controlled breath type enables the patient to take spontaneous breaths, and pressure is automatically adjusted to encourage the patient to pull the targeted tidal volume.

Trending Software: Shows time-stamped displays of ventilator settings up to 53 parameters of monitored patient data and specific events for intervals up to 72 hours.

Respiratory Mechanics Software: Enables monitoring of key respiratory parameters for easy assessment of patient status.

Tube Compensation Software: Accurately overcomes the work of breathing imposed by the artificial airway. Helps clinicians better recognize when a patient is ready to extubate.

Features Designed for Safety

- **Circuit disconnect detection:** Monitors circuit pressure and effective patient volume to promptly detect circuit disconnection.
- **Automatic patient detection:** Helps prevent inadvertent changes from existing ventilator settings to standby modes that do not provide ventilation.
- **Ongoing background checks:** Assesses the proper function of the ventilator's electronics and pneumatics hardware continuously during ventilation.
- **Heated expiratory filters:** Traps 99.97% of pathogens to shield patients, clinicians and visitors from exposure to viruses and bacteria from exhaled gases.^{18,19}
- **Puritan Bennett™ 803 extended backup power source (BPS):** Minimizes risks during power outages by providing up to four hours of continuous battery backup power.²⁰



Accessories

- **Puritan Bennett™ 840 ventilator pole cart:** Enables the ventilator to fit into confined spaces when a compressor is not required. Features dual-wheel castors and an ergonomic 360-degree handle system to improve maneuverability.
- **Puritan Bennett™ 840 compressor mount cart:** Offers a sturdy compressor mount plus either a one-hour or four-hour backup power supply.
- **Puritan Bennett™ 806 compressor:** Fits the compressor mount cart to provide a safe alternative air source if bottled or wall air are not available.
- **Puritan Bennett™ 840 ventilator cart:** Offers lightweight, easy maneuverability in a cost-effective package.



ORDERING INFORMATION

Standard Accessories		Catalog Number
Flex arm		4-032006-00
Inspiratory bacteria filter		
Disposable filter (D/Flex, carton of 12)		4-074601-00
Expiratory bacteria filter and collector vial		
Disposable filter (D/X800, carton of 12)		4-076887-00
Test hose		4-018506-00
Test lung		4-000612-00
Oxygen hose assembly, DIS (U.S.)		4-001474-00
Air hose assembly, DIS (U.S.)		4-006541-00
Power cord (North America)		4-071420-00
Operator's and technical reference manual		
English		4-075609-00
Spanish		4-070147-00
Software Options		
NeoMode 2.0 Software Option		10051492
Leak Compensation Software Option		10035870
Tube Compensation Software Option		4-076371-00
Bi-Level Software Option		4-076064-00
Volume Ventilation Plus Software Option		4-078126-00
PAV™*+ Option Kit		4-078203-00
Respiratory Mechanics Option Kit		10019218
Trending Software Option		10020408
Optional Accessories		
Service manual, English		4-070089-00
Puritan Bennett™ 840 Ventilator Cart and Accessories		
Puritan Bennett™ 840 Ventilator Cart with 1 Hr BPS		10000193
Puritan Bennett™ 840 Ventilator Cart with 4 Hr BPS		10000194
Wall-Air Water Trap Kit		4-075315-00
Fisher & Paykel Humidifier Mounting Kit for Black Carts		4-075313-00
Puritan Bennett™ 803 Backup Power Supply (BPS)		10030274
Battery Replacement Kit		4-070523-SP
Puritan Bennett™ 840 Ventilator Compressor Mount Cart and Accessories		
Puritan Bennett™ 840 Ventilator Compressor Mount Cart with 1 Hr BPS		10046822
Puritan Bennett™ 840 Ventilator Compressor Mount Cart with 4 Hr BPS		10046823
Cylinder Mount Bracket Kit		10045586
Wall-Air Water Trap Kit		10045588
Universal Humidifier Mounting Bracket for White Carts		10045589

Puritan Bennett™ 840 Ventilator Pole Cart and Accessories	
Puritan Bennett™ 840 Ventilator Pole Cart and 1 Hr BPS	10046826
Puritan Bennett™ 840 Ventilator Pole Cart and 4 Hr BPS	10046827
Cylinder Mount Bracket Kit	10045578
Wall-Air Water Trap Kit	10045588
Humidifiers and Breathing Circuits	
Reusable, adult, with heated wire, for Fisher & Paykel	G-061235-00
Reusable, adult, without heated wire	G-061208-SP
Reusable, pediatric, with heated wire, for Fisher & Paykel	G-061237-00
Reusable, pediatric, without heated wire	G-061223-00
Humidifiers and Breathing Circuits (continued)	
Humidifier Base	4-MR850-00
Puritan Bennett™ MR 850 Starter Kit	4-070773-00
Fisher & Paykel Mount Kit for Universal Bracket	10081874
Hudson Mount Kit for Universal Bracket	10081875
Kendall A2000 Mount for Universal Bracket	10081876
Fisher & Paykel Shroud Kit	10081877
Mounting Plate - Fisher & Paykel Shroud	10081785
Inspiratory bacteria filter	
Reusable filter (Re/Flex, each)	4-074600-00
Neo filter and adapter	
Neo disposable filter (carton of 12)	4-076408-00
Neo filter adapter	4-076405-00
Expiratory bacteria filter and collector vial	
Reusable filter (Re/X800, each)	4-070305-00
Reusable collector vial (Re/X800, each)	4-074647-00
Drain Bag and Drain Bag Accessories	
Drain bag, disposable (package of 25)	4-048491-00
Drain bag tubing, disposable (package of 10)	4-048493-00
Clamp, reusable (package of 5)	4-048492-00
Drain cap	4-074613-00
Filters and Sensors	
Seal, expiratory filter	4-070311-00
Filter, foam, compressor inlet	4-074374-00
Oxygen sensor††	4-072214-00
Nebulizer	
Aeroneb®* Pro Nebulizer	4-AP6000-US
Preventive Maintenance Kits	
10,000-hour preventive maintenance kit††, BDU/GUI	4-078179-00
15,000-hour preventive maintenance kit††, compressor	4-076805-00

†† Oxygen sensor to be replaced every two years or as necessary by a qualified service technician. Preventive maintenance kits must be installed by a qualified service technician.



References

1. Siegel MD. Management of agitation in the intensive care unit. *Clin Chest Med.* 2003;24(4):713-725.
2. Tate JA, Devito Dabbs A, Hoffman LA, Milbrandt E, Happ MB. Anxiety and agitation in mechanically ventilated patients. *Qual Health Res.* 2012;22(2):157-173.
3. Thille AW, Rodriguez P, Cabello B, Lellouche F, Brochard L. Patient-ventilator asynchrony during assisted mechanical ventilation. *Intensive Care Med.* 2006;32(10):1515-1522.
4. de Wit M, Miller KB, Green DA, Ostman HE, Gennings C, Epstein SK. Ineffective triggering predicts increased duration of mechanical ventilation. *Crit Care Med.* 2009;37(10):2740-2745.
5. Epstein SK. Optimizing patient-ventilator synchrony. *Semin Respir Crit Care Med.* 2001;22(2):137-152.
6. Levine S, Nguyen T, Taylor N, et al. Rapid disuse atrophy of diaphragm fibers in mechanically ventilated humans. *N Engl J Med.* 2008;358(13):1327-1335.
7. Hermans G. Increased duration of mechanical ventilation is associated with decreased diaphragmatic force: a prospective observational study. *Crit Care.* 2010;14:R127.
8. Pohlman MC, et al. Excessive tidal volume from breath stacking during lung-protective ventilation for acute lung injury. *Crit Care Med.* 2008;36(11):3019-3023.
9. Anzueto A, Peters JJ, Tobin MJ, et al. Effects of prolonged controlled mechanical ventilation on diaphragmatic function in healthy adult baboons. *Crit Care Med.* 1997;25(7):1187-1190.
10. Wilkins RL, Stoller JK, Scanlan CL. *Egan's Fundamentals of Respiratory Care.* 8th ed. Louis, MO: Mosby; 2003.
11. Haitsma JJ. Diaphragmatic dysfunction in mechanical ventilation. *Curr Opin Anaesthesiol.* 2011;24(2):214-218.
12. Puritan Bennett™ 840 ventilator operations manual
13. Costa R, Spinazzola G, Cipriani F, et al. A physiologic comparison of proportional assist ventilation with load-adjustable gain factors (PAV+) versus pressure support ventilation (PSV). *Intensive Care Med.* 2011;37(9):1494-1500.
14. Xirouchaki N, Kondili E, Klimathianaki M, Georgopoulos D. Is proportional-assist ventilation with load-adjustable gain factors a user-friendly mode? *Intensive Care Med.* 2009;35(9):1599-1603.
15. Xirouchaki N, Kondili E, Vaporidi K, et al. Proportional assist ventilation with load-adjustable gain factors in critically ill patients: comparison with pressure support. *Intensive Care Med.* 2008;34(11):2026-2034.
16. Bosma K, Ferreyra G, Ambrogio C, et al. Patient-ventilator interaction and sleep in mechanically ventilated patients: pressure support versus proportional assist ventilation. *Crit Care Med.* 2007;35(4):1048-1054.
17. Younes M. Proportional assist ventilation, a new approach to ventilatory support. *Theory. Am Rev Respir Dis.* 1992;145(1):114-120.
18. Younes M, et al. Proportional Assist Ventilation. In: Tobin M. *Principles and Practice of Mechanical Ventilation.* McGraw-Hill. 2006: 335-364.
19. Younes M, Webster K, Kun J, Roberts D, Masiowski B. A method for measuring passive elastance during proportional assist ventilation. *Am J Respir Crit Care Med.* 2001;164(1):50-60.
20. Grasso S, Ranieri WM, Brochard L, et al. Closed loop proportional assist ventilation (PAV): Results of a phase II multicenter trial. *Am J Respir Crit Care Med.* 2001, 163:A303.
21. Younes M, Riddle W, Polacheck J. A model for the relationship between respiratory neural and mechanical outputs: III. Validation. *J Appl Physiol.* 1981;51(4):990-1001.



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG.™* and ®* are trademarks of their respective owner. Proportional Assist and PAV are registered trademarks of The University of Manitoba, Canada. Used under license. Other brands are trademarks of a Covidien company. © 2013 Covidien. 13-VE-0013 VE27010

6135 GUNBARREL AVENUE
BOULDER, CO
80301

800-635-5267

WWW.COVIDIEN.COM