

Great. What we need the following: I appreciate your help.

1. tube scans on the tube 90,025
2. what kind of tube model number and age CXB-750D/4A
3. mas on the tube. 60KV, 150KV, 400mA, 7.5 million H.U.
4. Options.

- ① Voltage Regulator (URDU)
- ② Licom MWM (COT-32D)
- ③ Recon Unit, Real-Time (CCFR-007A.00A)
- ④ Display Console PC (FRUX74-000211.023)
- ⑤ PCP Study Split (COT-44A/3B)